

#### Fluorochromes:

- **Propidium iodide (PI):** Prepare a stock solution of 1 mg/ml and filter through a 0.22 µm filter. Store in 1 ml aliquots at -20°C. The working concentration of PI is usually 50 µg/ml.
- **4',6-diamidino-2-phenylindole (DAPI):** Prepare stock solution of 0.1 mg/ml and filter through a 0.22 µm filter. Store in 1 ml aliquots at -20°C. The working concentration of DAPI is usually 4 µg/ml.

#### Isolation buffers:

Isolation buffers must be prepared using either single- or double-distilled water, filtered through a 0.22 µm filter to remove suspended particles, and stored as specified. The pH of the buffers is adjusted either with 1M NaOH or with 1N HCl.

- **Lysis buffer LB01 (Doležel *et al.*, 1989):** 15 mM Tris, 2 mM Na<sub>2</sub>EDTA, 0.5 mM spermine.4HCl, 80 mM KCl, 20 mM NaCl, 0.1 % (v/v) Triton X-100. Adjust to pH 7.5. Add β-mercaptoethanol to give a final concentration of 15 mM. Store the buffer either at 4°C if used regularly or at -20°C in 10 ml aliquots. Concentrations of 0.5%, 4% and 8% of Triton X-100 can be used in recalcitrant materials.
- **Tris MgCl<sub>2</sub> buffer (Pfosser *et al.*, 1995):** 200 mM Tris, 4 mM MgCl<sub>2</sub>, 0.5 % (v/v) Triton X-100. Adjust pH to 7.5 and store at 4°C.
- **Galbraith's buffer (Galbraith *et al.*, 1983):** 45 mM MgCl<sub>2</sub>, 20 mM MOPS, 30 mM sodium citrate, 0.1 % (v/v) Triton X-100. Adjust pH to 7.0. Store the buffer either at 4°C if used regularly or at -20°C in 10 ml aliquots.
- **General purpose buffer (Loureiro *et al.*, 2007):** 0.5 mM spermine.4HCl, 30 mM sodium citrate, 20 mM MOPS, 80 mM KCl, 20 mM NaCl, 0.5 % (v/v) Triton X-100. Adjust to pH 7.0. Store the buffer either at 4°C if used regularly or at -20°C in 10 ml aliquots.

## Checklist for plant genome size estimation by flow cytometry GSAD: a genome size in the Asteraceae database

- **Woody plant buffer (Loureiro *et al.*, 2007):** 200 mM Tris, 4 mM MgCl<sub>2</sub>, 2 mM Na<sub>2</sub>EDTA, 86 mM NaCl, 10 mM sodium metabisulphite, 1 % PVP-10, 1 % (v/v) Triton X-100. Adjust to pH 7.5. Store the buffer either at 4°C if used regularly or at -20°C in 10 ml aliquots.
- **MgSO<sub>4</sub> buffer (Arumuganathan and Earle, 1991):** 9.53 mM MgSO<sub>4</sub>, 47.67 mM KCl, 4.77 mM HEPES, 6.48 mM DTT (*see note 1*), 0.25 % (v/v) Triton X-100. Adjust to pH 8.0. Store the buffer either at 4°C if used regularly or at -20°C in 10 ml aliquots.
- **De Laat's buffer (de Laat and Blaas, 1984):** 15 mM HEPES, 1mM Na<sub>2</sub>EDTA, 0.2 % (v/v) Triton X-100, 80 mM KCl, 20 mM NaCl, 15 mM DTT, 0.5 mM spermine.4HCl, 300 mM sucrose. Adjust to pH 7.0 and store at 4°C.
- **Ebihara's buffer (Ebihara *et al.*, 2005):** 50 mM Na<sub>2</sub>SO<sub>3</sub>, 50 mM Tris, 40 mg/mL PVP-40 (*see Note 7*), 140 mM β-mercaptoethanol. Adjust to pH 7.5 and store at 4°C.
- **Seed buffer (Matzk *et al.*, 2001):** 5 mM MgCl<sub>2</sub>, 85 mM NaCl, 100 mM Tris, 0.1 % Triton X-100. Adjust to pH 7.0 and store at 4°C.
- **Otto's buffers (Otto, 1992):**  
Otto I: 100 mM citric acid monohydrate, 0.5 % (v/v) Tween 20 (cell culture tested grade of Tween 20 from Sigma-Aldrich (cat. no. P2287) is used. Tween 20 for molecular biology (Sigma cat. no. P9416) is not suitable for FCM. Store at 4°C.  
  
Otto II: 400 mM Na<sub>2</sub>HPO<sub>4</sub>. Store at room temperature. The fluorochrome (DAPI or PI, *see above*) can be added to Otto II before adjusting the final volume of the stock solution. If this is done the buffer should be stored in the dark at room temperature.
- **Baranyi's buffer (Baranyi and Greilhuber, 1995):** Baranyi's solution I: 100 mM citric acid monohydrate, 0.5 % (v/v) Triton X-100. Store at 4°C.  
  
Baranyi's solution II: 400 mM Na<sub>2</sub>HPO<sub>4</sub>, 10 mM sodium citrate, 25 mM sodium sulphate. Store at room temperature.

## Checklist for plant genome size estimation by flow cytometry GSAD: a genome size in the Asteraceae database

- **Mishiba's buffer (Mishiba *et al.*, 2000):** Solution A: (see recipe for Galbraith buffer, i.e. buffer 3 above). Solution B: 10 mM Tris, 50 mM sodium citrate, 2 mM MgCl<sub>2</sub>, 1 % PVP-40 (original recipe used PVP K-30), 0.1 % Triton X-100, 18 mM β-mercaptoethanol. Adjust to pH 7.5. Store at 4°C.
- **Brown's Nuclear Buffer:** 45 mM MgCl<sub>2</sub>, 60 mM sodium citrate, 20 mM 4-morpholinepropane sulfonate pH7, 0.1% (w/v) Triton X-100, 1% polyvinyl pyrrolidone (~10,000 M<sub>r</sub>, Sigma P6755), 5 mM sodium metabisulphite (syn. pyrosulphite, Na<sub>2</sub>S<sub>2</sub>O<sub>5</sub> M<sub>r</sub> 190; added every twelve hours from 1 M frozen stock). This is modified from Galbraith *et al.* (1983); stronger buffer (for acidic plants or fruit), protectant against tanning (PVP) and antioxidant (metabisulphite, not toxic like β-mercaptoethanol) are used. Triton may be raised to 0.5% for oily tissues ( *Pistacia* sp. etc.) and to lyse chloroplasts which fluoresce red; keep 10% (w/v) Triton X-100 stock, autoclaved and stored at 4°C for this purpose.
- **Doležel's Nuclear Buffer (Doležel *et al.*, 2007):** 20 mM NaCl, 80 mM KCl, 20 mM MgSO<sub>4</sub>, 2 mM EDTA.Na<sub>2</sub>, 0.5 mM spermine.HCl, 15 mM Tris pH 7.5, 0.1% (w/v) Triton X-100, 15 mM β-mercaptoethanol (1 μl / ml) added daily – but this is toxic. It is better to use metabisulphite.

### References

- Arumuganathan K, Earle ED. 1991. Nuclear DNA content of some important plant species. *Plant Molecular Biology Reporter*, 9: 208–218.
- Baranyi M, Greilhuber J. 1995. Flow cytometric analysis of genome size variation in cultivated and wild *Pisum sativum* (Fabaceae). *Plant Systematics and Evolution*, 194: 231–239.
- de Laat AMM, Blaas J. 1984. Flow cytometric characterization and sorting of plant chromosomes. *Theoretical and Applied Genetics*, 67: 463–467.
- Doležel J, Binarova P, Lucretti S. 1989. Analysis of nuclear DNA content in plant cells by flow cytometry. *Biologia Plantarum*, 31: 113–120.



## Checklist for plant genome size estimation by flow cytometry GSAD: a genome size in the Asteraceae database

Doležel J, Greilhuber J, Suda J. 2007. Flow Cytometry with Plant Cells Analysis of Genes, Chromosomes and Genomes. Weinheim: Wiley-VCH Verlag GmbH & Co. KGaA.

Ebihara A, Ishikawa H, Matsumoto S *et al.* 2005. Nuclear DNA, chloroplast DNA, and ploidy analysis clarified biological complexity of the *Vandenboschia radicans* complex (Hymenophyllaceae) in Japan and adjacent areas. *American Journal of Botany*, 92: 1535–1547.

Galbraith DW, Harkins KR, Maddox JM, Ayres NM, Sharma DP, Firoozabady E. 1983. Rapid flow cytometric analysis of the cell cycle in intact plant tissues. *Science*, 220: 1049–1051.

Loureiro J, Rodriguez E, Doležel J, Santos C. 2007. Two new nuclear isolation buffers for plant DNA flow cytometry: a test with 37 species. *Annals of Botany*, 100: 875–888.

Matzk F, Meister A, Brutovská R, Schubert I. 2001. Reconstruction of reproductive diversity in *Hypericum perforatum* L. opens novel strategies to manage apomixis. *The Plant Journal*, 26: 275–282.

Mishiba KI, Ando T, Mii M *et al.* 2000. Nuclear DNA content as an index character discriminating taxa in the genus *Petunia* sensu Jussieu (Solanaceae). *Annals of Botany*, 85: 665–673.

Otto F. 1992. Preparation and staining of cells for high-resolution DNA analysis. In: Radbruch A ed. Flow cytometry and cell sorting. Berlin: Springer-Verlag, 101–104.

Pfossen M, Amon A, Lelley T, Heberlebers E. 1995. Evaluation of sensitivity of flow-cytometry in detecting aneuploidy in wheat using disomic and ditelosomic wheat-rye addition lines. *Cytometry*, 21: 387–393.