

	FISH PROTOCOL	°C	Time	Machinery / tools
FIRST DAY - MORNING	<i>Go and find ice. Set incubator at 37°C with a humid chamber (wet paper with 2xSSC) inside it. Make sure that there are enough plastic (i.e. parafilm) coverslips.</i>			
PREVIOUS (facultative)	Rinse slides with EtOH 95% 10 min, air-dry.	RT		Hellendhal jars
RNAse TREATMENT (100 µg/ml)	For 4 slides: 8 µl Rnase (10 mg/ml) + 792 µl 2xSSC . Vortex.	Ice		Vortex Pipette
	Incubate slides with 200 µl RNAse with plastics (i.e. parafilm) in humid chamber.	37	1 h	Incubator Humid chamber Vortex Pipette
WASHING	In 2xSSC with gentle agitation, remove plastic coverslips with forceps (or allow the to fall off).	RT	5 min x 3	3D-Shaker Hellendhal jars
Prior to pepsin TREATMENT	Treat slides with HCl 0.01N .	RT	2 min	Hellendhal jars Pipette
	Dilute the pepsin stock (10 mg/ml) 1:200 in 0.01N HCl.	Ice		Centrifuge Pipette
Pepsin TREATMENT (0.05 mg/ml)	Incubate slides with 100 µl pepsin 0.05 mg/ml , with plastic coverslips (in h. ch.)	37	15 min	Pipette
WASHING	Remove plastic coverslips. Wash with ultrapure water (sterilised).	RT	2 min	Hellendhal jars
	2xSSC	RT	5 min x 2	3D-Shaker Hellendhal jars
FIXATION (facultative depending on the materials, but usually recommended)	Fix slides in neutral formaldehyde 3.7% .	RT	10 min	Hood
DEHYDRATION	In ethanol series: EtOH at 50%, 70% i 100% (at -20°C).	RT	2 min each	3D-Shaker Hellendhal jars
DRYING	Air drying.	RT	>1h 30min	

Everything in red is either optional or alternative.

	FISH PROTOCOL	°C	Time	Machinery / utiles
FIRST DAY - AFTERNOON	<i>Go and find ice. Prepare baths, one at 75°C and one at 55°C. Preheat humid chamber at 37°C.</i>			
PROBE MIXTURE PREPARATION	For 30 µl per slide mix (x4 slides): -formamide: 15 µl (60 µl) -50% dextran sulfate: 6 µl (24 µl) -20xSSC: 3 µl (12 µl) -labelled rDNA probe 1 : 50 to 100 ng per slide -labelled rDNA probe 2 : 50 to 100 ng per slide - water: up to 30 µl (up to 120 µl)	Ice		Eppendorfs Hood Pipette
PROBE MIXTURE DENATURATION	Denature probe mixture.	75	10 min	Bath 75
	Cool on ice bath, rapidly. Vortex and spin.	Ice		Ice Vortex Centrifuge
	Incubate probe mixture. Spin again.	Ice	5 min	Centrifuge
BEGINNING HYBRIDISATION PROCESS	Apply 30 µl per slide of the denatured probe per slide. Cover with plastic coverslips. Put slides into closable and <u>preheated</u> boxes (at 75°C) with <u>thin bottom</u> (cleant previously with ethanol), and these in the bath.	75	10 min	Bath 75 Box for bath Pipette
CONTINUING HYBRIDISATION PROCESS	Transfer boxes CAREFULLY to the other bath (at 55°C).	55	5 min	Bath 55 Boxes for bath
HYBRIDISATION PROCESS	Transfer slides CAREFULLY to a <u>preheated</u> (37°C) humid chamber (wet paper with 2xSSC). Seal the chamber with parafilm.	37	>18 h	Incubator Humid chamber Parafilm

Everything in red is either optional or alternative.

	FISH PROTOCOL	°C	Time	Machinery / tools
SECOND DAY - MORNING				
WASHING	2xSSC	42	5 min x 2	Bath with linear shaking Hellendhal jars
ASTRINGENT WASHING	0.1xSSC	42	5 min x 2	Bath with linear shaking Hellendhal jars
ALTERNATIVE ASTRINGENT WASHING	50% formamide in 1xSSC	42	10 min	Bath with linear shaking Hellendhal jars
WASHING	2xSSC	42 cool to 20	5 min x 2	Bath with linear shaking Hellendhal jars
WASHING	2xSSC	RT	5 min	Hellendhal jars
WASHING	4xSSCT	RT	7 min	Hellendhal jars
RINSE	Briefly rinse in 1xPBS and dry	RT		Hellendhal jars
DRYING	Allow to dry briefly (not too wet, not too dry)	RT		
MOUNTING	Mounting in Vectashield (12 µl per slide) . Allow Vectashield to expand (in the dark).	RT, dark	5-10 min	Big coverslips (24 x 40 mm) Pipette
STORAGE	In appropriate folders (in the dark).	4	for ever	Fridge Folders

Everything in red is either optional or alternative.